

## Quality Improvements of Xe<sup>+</sup> PFIB Prepared Surfaces of Aluminium

Van Leer, B.<sup>1</sup> and Passey, R.<sup>1</sup>

<sup>1</sup> Thermo Fisher Scientific, United States

Advanced engineering of materials containing aluminium requires a critical understanding of the material's microstructure. Cross-section and S/TEM sample preparation of aluminium and aluminium alloys to characterize grain boundaries or prepare surfaces for analytical work such as EBSD by focused ion beam (FIB) continues to be a major interest in metallurgical analysis because of FIB's ability to quickly prepare site specific specimens while eliminating damage from mechanical polishing or electro-polishing [1]. However, it is well known that gallium is particularly mobile in aluminium and its alloys, often diffusing along grain boundaries and interfaces, which can lead to embrittlement [2,3]. One attempt to mitigate this artifact is to subsequently polish the Ga<sup>+</sup> FIB prepared specimens in a board beam polisher to improve the sample's surface quality.


Characterization of large grain aluminium samples by Ga<sup>+</sup> FIB requires longer processing time when investigating multiple grains. Recent instrumentation using plasma FIB (PFIB) technology and Xe<sup>+</sup> ions offers an alternative ion species for FIB milling with increased milling rates because of its ability to deliver 30 - 40 times more current compared to Ga<sup>+</sup> FIBs. While the measured sputter rate of aluminum using Ga<sup>+</sup> and Xe<sup>+</sup> differs by about 25% (0.31  $\mu\text{m}^3/\text{nC}$  [Ga] and 0.41  $\mu\text{m}^3/\text{nC}$  [Xe]), the ability to prepare gallium free thin sections for S/TEM analysis or large area cross-sections offers a solution to FIB milling artifacts seen with a Ga<sup>+</sup> FIB. This also improves statistical analysis of grain populations.

To compare and contrast the effects of xenon ions to gallium ions, FIB prepared cross-sections and thin S/TEM specimens of commercial grade 6061 T6 aluminum were prepared using the Helios G4 DualBeam™ using 30 kV Ga<sup>+</sup> ions and a Helios G4 PFIB DualBeam using 30 kV Xe<sup>+</sup> ions. Sample quality of the cross-sections were evaluated by comparing EBSD index rates and band contrast. Sample quality of the S/TEM specimens were evaluated by measuring FIB sidewall damage and observation of Ga contamination at grain boundaries in a Themis Z™ TEM operating at 300 keV.


A comparison of FIB prepared cross-sections (30 kV; 15 nA) in 6061 T6 reveals a 25% increase in surface quality when using xenon ions. Figure 1 shows inverse pole figure maps overlaid with SEM images. A comparison of FIB sidewall damage between specimens prepared in Ga<sup>+</sup> FIB and Xe<sup>+</sup> PFIB revealed a decrease in sidewall damage by as much as 30% depending on the accelerating voltage applied. Chart 1 reveals FIB sidewall damage in 6061 T6 for xenon and gallium at 30 kV, 5 kV and 2 kV respectively.

### References:

- [1] L.A. Giannuzzi et al., Mater. Res. Soc. Symp. Proc. **480** (1997) p 19.
- [2] K.A. Unocic et al., Journal of Microscopy **240** (2010) p 227.
- [3] R.C. Hugo and R.G. Hoagland, Scr. Mater. **41** (1999) p 1341.



**Figure 1.** Inverse pole figure results of aluminium 6061 T6 with Ga<sup>+</sup> prepared surface (left) and Xe<sup>+</sup> prepared surface (right).



**Chart 1.** FIB sidewall damage in aluminium 6061 T6 at various accelerating voltages.